NKMAXBIO We support you, we believe in your research

Human CNDP2/CPGL antibody

Catalog Number: ATGA0562

PRODUCT INFORMATION

Catalog number

ATGA0562

Clone No.

AT15E5

Product type

Monoclonal antibody

UnitProt No.

Q96KP4

NCBI Accession No.

NP 060705

Alternative Names

Cytosolic non-specific dipeptidase, CN2, CPGL, HsT2298, PEPA

Additional Information

This product was produced from tissue culture supernatant.

PRODUCT SPECIFICATION

Antibody Host

Mouse

Reacts With

Human

Concentration

1mg/ml (determined by BCA assay)

Formulation

Liquid in. Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% glycerol

Immunogen

Recombinant human CNDP2 (1-475aa) purified from E. coli

Isotype

IgG1 kappa

Purification Note

By protein-A affinity chromatography

Application

ELISA, WB, ICC/IF

Usage

The antibody has been tested by ELISA, Western blot and ICC/IF analysis to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.



NKMAXBIO We support you, we believe in your research

Human CNDP2/CPGL antibody

Catalog Number: ATGA0562

Storage

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description


CNDP2, known as tissue carnosinase and peptidase A, is a nonspecific dipeptidase rather than a selective carnosinase and is a secreted peptidase homologous to M20 peptidases. Also known as glutamate carboxypeptidase-like protein 1 (CPGL), CNDP2 is ubiquitously expressed throughout human tissue. CNDP2 Hydrolyzes a variety of dipeptides including L-carnosine but has a strong preference for Cys-Gly and may be play a role as tumor suppressor in hepatocellular carcinoma (HCC) cells.

General References

Zhang, P., et al. (2006) Cancer res 12: 6617-6625. Kaur, H., et al. (2009) J.Biol.Chem 284: 14493-14502. Licker, V., et al. (2012) J proteomics 375(15): 4656-67.

DATA

Western blot analysis (WB)



The tissue lysate (40ug) was resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human CNDP2/CPGL antibody (1:1000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.

Lane 1.: Mouse kidney tissue lysate


Immunocytochemistry/Immunofluorescence (ICC/IF)



NKMAXBio We support you, we believe in your research

Human CNDP2/CPGL antibody

Catalog Number: ATGA0562



ICC/IF analysis of CNDP2/CPGL in HeLa cells. The cell was stained with ATGA0562 (1:100). The secondary antibody (green) was used Alexa Fluor 488. DAPI was stained the cell nucleus (blue).

