NKMAXBio We support you, we believe in your research

A431 Lysate

Catalog Number: ATGL0007

PRODUCT INFORMATION

Catalog number

ATGL0007

PRODUCT SPECIFICATION

Concentration

2mg/ml (determined by BCA assay)

Formulation

Liquid in. SDS-PAGE Sample Buffer (0.12M Tris-HCl buffer pH6.8, 0.4% SDS, 5% glycerol, 1% beta-mercaptoethanol, 0.02% Bromophenol blue)

Storage

Aliquot the lysate into single tube and then storing them -80C or below. Avoid repeated freezing and thawing cycles.

Preparation before usage

Heat lysate to 95C for 5 minutes and rapidly cooling before use.

General description


A431 cells were established from an epidermoid carcinoma in the skin/epidermis of an 85 year old female patient. A431 cells are used in studies of the cell cycle and cancer-associated cell signalling pathways since they express abnormally high levels of the EGFR. As such they are often used as a positive control for EGFR expression.

Quality Control

Western Blot analyzed by both GAPDH and Beta tubulin antibody

DATA

Western blot analysis (WB)



All lanes: A431 Whole Cell Lysate (40ug)

Lane 1. : Anti-GAPDH antibody (ATGA0181) Lane 2. : Anti-Beta tubulin antibody (ATGA0196)

