NKMAXBIO We support you, we believe in your research

Recombinant human S100B protein

Catalog Number: ATGP0628

PRODUCT INFORMATION

Expression system

E.coli

Domain

1-92aa

UniProt No.

P04271

NCBI Accession No.

NP 006263

Alternative Names

S100 calcium binding protein B, S100 calcium binding protein beta neural, S100beta, S100 protein beta chain, S100 protein subunit beta, S100 calcium-binding protein B

PRODUCT SPECIFICATION

Molecular Weight

12.8 kDa (112aa) confirmed by MALDI-TOF

Concentration

1mg/ml (determined by Bradford assay)

Formulation

Liquid in. 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol

Purity

> 90% by SDS-PAGE

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

S100B is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. This protein function in Neurite extension, proliferation of melanoma cells, stimulation of Ca2+ fluxes, inhibition of PKC-mediated phosphorylation, astrocytosis and axonal proliferation, and inhibition of microtubule assembly. Recombinant human S100B



NKMAXBio We support you, we believe in your research

Recombinant human S100B protein

Catalog Number: ATGP0628

protein, fused to His-tag at N-terminus, was expressed in E. coli and purified by using conventional chromatography techniques.

Amino acid Sequence


MGSSHHHHHH SSGLVPRGSH MSELEKAMVA LIDVFHQYSG REGDKHKLKK SELKELINNE LSHFLEEIKE QEVVDKVMET LDNDGDGECD FQEFMAFVAM VTTACHEFFE HE

General References

van Dieck J., et al. (2010) Oncogene. 29(14):2024-35. Ayilavarapu S., et al. (2010) J Immunol. 184(3):1507-15.

DATA

SDS-PAGE



15% SDS-PAGE (3ug)

3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

