

Recombinant human GLIPR2 protein

Catalog Number: ATGP0689

PRODUCT INFORMATION

Expression system

E.coli

Domain

1-154aa

UniProt No.

Q9H4G4

NCBI Accession No.

NP_071738

Alternative Names

Golgi-associated plant pathogenesis-related protein 1, C9orf19, GAPR-1, Golgi-associated plant pathogenesis-related protein 1

PRODUCT SPECIFICATION

Molecular Weight

19.3 kDa (174aa) confirmed by MALDI-TOF

Concentration

1mg/ml (determined by Bradford assay)

Formulation

Liquid in. 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 1mM DTT, 0.1M NaCl.

Purity

> 90% by SDS-PAGE

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

GLIPR2 is closely related to plant pathogenesis-related (PR-1) proteins, which are upregulated in response to pathogen attack. This protein is found within lipid-enriched microdomains on the cytosolic side of the endomembrane system. GLIPR2 is tightly anchored to membranes and absent from the cytosol, although it does not possess a membrane-spanning domain. Recombinant human GLIPR2 protein, fused to His-tag at N-terminus, was expressed in E. coli and purified by using conventional chromatography techniques.

Recombinant human GLIPR2 protein

Catalog Number: ATGP0689

Amino acid Sequence


<MGSSHHHHHH SGLVPRGSH> MGKSASKQFH NEVLKAHNEY RQKHGVPPLK LCKNLNREAQ QYSEALASTR
ILKHSPSSSR GQCGENLAWA SYDQTGKEVA DRWYSEIKNY NFQQPGFTSG TGHFTAMVWK NTKKMGVGKA
SASDGSSFVV ARYFPAGNVV NEGFFEENVL PPKK

General References

Eberle HB., et al. (2002) J Cell Sci. 115(4):827-38.

DATA

SDS-PAGE



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain