NKMAXBIO We support you, we believe in your research

Recombinant rat IL-1 alpha/IL1A protein

Catalog Number: ATGP3192

PRODUCT INFORMATION

Expression system

E.coli

Domain

115-270aa

UniProt No.

P16598

NCBI Accession No.

NP 058715

Alternative Names

Interleukin-1 alpha precursor, Interleukin-1 alpha, IL1F1, IL-1A, IL-1 alpha

PRODUCT SPECIFICATION

Molecular Weight

20.2 kDa (179aa) confirmed by MALDI-TOF

Concentration

1mg/ml (determined by absorbance at 280nm)

Formulation

Liquid in. Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Purity

> 95% by SDS-PAGE

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

IL1A also known as interleukin-1 alpha. IL1A belongs to the interleukin 1 cytokine family. IL1A is part of Interleukin-1. IL1 stimulates thymocyte proliferation by inducing IL2 release, B-cell maturation and proliferation, and fibroblast growth factor activity. IL1 proteins are involved in the inflammatory response, especially IL1A is an early and important mediator of inflammation after injury. Recombinant rat IL1A, fused to His-tag at N-terminus was expressed in E. coli and purified by using conventional chromatography techniques.



NKMAXBio We support you, we believe in your research

Recombinant rat IL-1 alpha/IL1A protein

Catalog Number: ATGP3192


Amino acid Sequence

MGSSHHHHHH SSGLVPRGSH MGSSAPHSFQ NNLRYKLIRI VKQEFIMNDS LNQNIYVDMD RIHLKAASLN DLQLEVKFDM YAYSSGGDDS KYPVTLKVSN TQLFVSAQGE DKPVLLKEIP ETPKLITGSE TDLIFFWEKI NSKNYFTSAA FPELLIATKE QSQVHLARGL PSMIDFQIS

General References

Brough D., et al. (2015) IUBMB Life. doi: 10.1002 Hermann E., et al. (1985) Z Rheumatol. 44(5): 207-12.

DATA



15% SDS-PAGE (3ug)

3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

