# NKMAXBIO We support you, we believe in your research

# **Human beta-Actin antibody**

Catalog Number: ATGA0606

### **PRODUCT INFORMATION**

# Catalog number

ATGA0606

#### Clone No.

AT4G6

### **Product type**

Monoclonal Antibody

#### UnitProt No.

P60709

#### **NCBI Accession No.**

NP 001092

#### **Alternative Names**

Actin, cytoplasmic 1, Beta-actin

#### **Additional Information**

This product was produced from tissue culture supernatant.

# **PRODUCT SPECIFICATION**

#### **Antibody Host**

Mouse

#### **Reacts With**

Human

# Concentration

0.2mg/ml (determined by BCA assay)

#### **Formulation**

Liquid in. Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% glycerol

### **Immunogen**

Recombinant human ACTB (1-375aa) purified from E. coli

### Isotype

IgG1 kappa

### **Purification Note**

By protein-A affinity chromatography

# **Application**

ELISA, WB, ICC/IF

#### **Usage**

The antibody has been tested by ELISA, Western blot and ICC/IF analysis to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.



# NKMAXBio We support you, we believe in your research

# **Human beta-Actin antibody**

Catalog Number: ATGA0606

#### **Storage**

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

#### **BACKGROUND**

#### **Description**


Beta-actin (human gene and protein symbol ACTB/ACTB) is one of six different actin isoforms which have been identified in humans. This is one of the two nonmuscle cytoskeletal actins. Actins are highly conserved proteins that are involved in cell motility, structure and integrity. Alpha actins are a major constituent of the contractile apparatus.


#### **General References**

Gunning PW., et al. (2015) Journal of Cell Science. 128(11): 2009-2019. Hanukoglu I., et al. (1983) Journal of Molecular Biology. 163(4): 673-8.

#### **DATA**

### Western blot analysis (WB)





The cell lysates (40ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human ACTB antibody (1:1000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.

Lane 1.: HeLa cell lysate Lane 2.: Jurkat cell lysate Lane 3.: PC3 cell lysate Lane 4.: HepG2 cell lysate Lane 5.: MCF7 cell lysate Lane 6.: U87-MG cell lysate Lane 7.: NIH-3T3 cell lysate

The mouse tissue lysates (40ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human ACTB antibody (1:1000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.

Lane 1.: Kidney tissue lysate Lane 2.: Spleen tissue lysate Lane 3.: Bone marrow tissue lysate Lane 4.: Stomach tissue lysate


Lane 5.: Liver tissue lysate



# NKMAXBio We support you, we believe in your research

# **Human beta-Actin antibody**


Catalog Number: ATGA0606



The HeLa cell lysates (40ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human ACTB antibody. Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.

Lane 1.: ACTB antibody 1:1000 Lane 2.: ACTB antibody 1:3000 Lane 3.: ACTB antibody 1:5000 Lane 4.: ACTB antibody 1:10000

# Immunocytochemistry/Immunofluorescence (ICC/IF)



ICC/IF analysis of ACTB in HeLa cells. The cell was stained with ATGA0570 (1:100). The secondary antibody (green) was used Alexa Fluor 488. DAPI was stained the cell nucleus (blue).

