

Anti-SARS-Cov-2 ELISA IgG

Catalog Number: ATGCF001

PRODUCT INFORMATION

Catalog number

ATGCF001

PRODUCT SPECIFICATION

Description

This Anti-SARS-CoV-2 ELISA (IgG) kit has been designed for the detection and the quantitative and qualitative determination of specific IgG antibodies against SARS-CoV-2 in plasma. This kit is intended for Research Use only and is not for use in diagnostic or therapeutic procedures. The Anti-SARS-CoV-2 ELISA (IgG) kit is an enzyme immunoassay (EIA). SARS-CoV-2 antigen is bound on the surface of the microtiter strips. Diluted patient plasma or standards are pipetted into the wells of the microtiter plate. A binding between the IgG antibodies of the plasma and the immobilized SARS-CoV-2 antigen takes place. After a one hour incubation at 37°C, the plate is rinsed with diluted wash solution, in order to remove unbound material. Then reconstituted conjugate is added and incubated for 1 hour at 37°C. After a further washing step, the substrate (TMB) solution is pipetted and incubated for 30 minutes. The color development is terminated by the addition of a stop solution, which changes the color from blue to yellow. The resulting dye is measured spectrophotometrically at the wavelength of 450nm. The concentration of the IgG antibodies is directly proportional to the intensity of the color.

[Kit insert_English version download link](#) [Kit insert_Korean version download link](#) [Brochure download link](#)

Storage

2-8°C

Product type

ELISA assay