# NKMAXBIO We support you, we believe in your research

## Recombinant human Fc gamma RIIB/CD32b protein

Catalog Number: ATGP4149

#### PRODUCT INFORMATION

## **Expression system**

**HEK293** 

#### **Domain**

46-217aa

#### UniProt No.

P31994

#### **NCBI Accession No.**

NP 001002274.1

#### **Alternative Names**

Low affinity immunoglobulin gamma Fc region receptor II-b isoform3, IgG Fc receptor II-b, CDw32, Fc-gamma RII-b, Fc-gamma-RIIb, FcRII-b, CD32, FCGR2B, FCG2, IGFR2, Fc fragment of IgG, low affinity II, low affinity IIb, Low affinity immunoglobulin gamma Fc region receptor II-c, IgG Fc receptor II-c, Fc-gamma RII-c, Fc fragment of IgG receptor IIb, Fc-gamma-RIIc, Fc gamma RIIb, CD32B, FcgammaRIIb, FCGR2, FCGR2C, FcGRIIB, FcRII-c

#### **Additional Information**

ATGP1822 has replaced with a catalog number ATGP4149.

## **PRODUCT SPECIFICATION**

#### **Molecular Weight**

20.1kDa (178aa)

#### Concentration

1mg/ml (determined by Absorbance at 280nm)

#### **Formulation**

Liquid in. Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

#### **Purity**

> 95% by SDS-PAGE

#### **Endotoxin level**

< 1 EU per 1ug of protein (determined by LAL method)

#### **Biological Activity**

Measured by its binding ability in a functional ELISA with Human IgG1 Fc (CAT# ATGP4094).

### Tag

His-Tag

## **Application**

SDS-PAGE, Bioactivity

## **Storage Condition**

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid



## Recombinant human Fc gamma RIIB/CD32b protein

Catalog Number: ATGP4149

repeated freezing and thawing cycles.

## **BACKGROUND**

#### **Description**

Fc gamma RIIB/CD32b, also known as FCGR2B, belongs to the Ig superfamily. It is a low-affinity receptor for the Fc region of immunoglobulins gamma that are complexed or aggregated. It participates in various effector and regulatory functions, such as phagocytosis of immune complexes and modulation of antibody production by B cells. Fc gamma RIIB/CD32b is expressed on B cells, monocytes, dendritic cells, neutrophils, mast cells, and basophils. When Fc gamma RIIB is ligated, it inhibits B cell expansion, plasma cell differentiation, production of autoimmune rheumatoid factors, and down-regulates TLR4 with reduced LPS responsiveness. Recombinant human Fc gamma RIIB/CD32b, fused to His-tag at C-terminus, was expressed in HEK293 cell and purified by using conventional chromatography techniques.

## **Amino acid Sequence**


APPKAVLKLE PQWINVLQED SVTLTCRGTH SPESDSIQWF HNGNLIPTHT QPSYRFKANN NDSGEYTCQT GQTSLSDPVH LTVLSEWLVL QTPHLEFQEG ETIVLRCHSW KDKPLVKVTF FQNGKSKKFS RSDPNFSIPQ ANHSHSGDYH CTGNIGYTLY SSKPVTITVQ AP<HHHHHHH>

#### **General References**

Van de Winkel, J. and P. Capes (1993) Immunol. Today 14:215-221. Raghaven, M. and P. Bjorkman (1996) Annu. Rev. Cell Dev. Biol. 12:181-220. Ravetch, J. and S. Bolland (2001) Annu. Rev. Immunol. 19:275-290.

## **DATA**

### SDS-PAGE




3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain

## **Biological Activity**

# Recombinant human Fc gamma RIIB/CD32b protein

Catalog Number: ATGP4149



Human Fc gamma RIIB/CD32b (CAT# ATGP4149) is coated at 10 ug/ml (100 ul/well) can bind Human IgG1 Fc (CAT# ATGP4094) in a Functional ELISA assay.

